
**LIGHTING &  
ELECTRONICS  
JENA**

# **LH-LED 100**

**LED-Lampenhous  
LED lamp housing**

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

Kollimator  
55 01 02-200-14


# **english**

The content of this document may change as design revisions of the product take place.  
It is then only the German version. The English part is not updated.

# LH-LED100

## LED lamp housing

OPERATING INSTRUCTIONS

- 
- © Distributing, duplicating in any form (by way of printing, copying, microfilming or by using any other technique) as well as processing this document by electronic means is permitted, unless explicitly approved. Violations will be prosecuted.

Subject to modifications in the interest of technical developments; these operating instructions are not subject to the revision service.

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Issue date: 06/26/2017  
Operating instructions No. 55 01 02-001-26A

## CONTENT

	Page
Copyright .....	24
Content .....	25
Illustrations and Tables .....	26
<b>1</b>	<b>Fundamentals .....27</b>
1.1	General instructions .....27
1.2	Instructions concerning the appliance's safety .....28
1.3	Shipping, storage and unpacking instructions .....30
1.4	Warranty instructions .....30
<b>2</b>	<b>Description of the appliance .....31</b>
2.1	Intended use .....31
2.2	Identification of the appliance .....31
2.3	Scope of supply .....31
2.4	Special features of the appliance .....31
<b>3</b>	<b>Technical data .....33</b>
<b>4</b>	<b>Putting the appliance into operation .....36</b>
4.1	Positioning the appliance .....36
4.2	Fastening the microscope adapter .....36
4.3	Connecting electrical equipment .....37
<b>5</b>	<b>Operation .....39</b>
5.1	Operating and display elements .....39
5.2	Collimator adjustment .....40
5.3	Setting the brightness level .....41
<b>6</b>	<b>Service .....42</b>
6.1	Maintenance and cleaning .....42
6.2	Spare parts & accessories .....42
6.3	Repairs .....42

---

## ILLUSTRATIONS AND TABLES

	Seite
Fig. 1: Luminous flux spectral distribution as compared with halogen lamps .....	34
Fig. 2: Fastening the microscope adapter .....	36
Fig. 3: Electrical connections .....	37
Fig. 4: View at the lamp housing's operating & display elements .....	39
Fig. 5: Collimator adjustment at the bottom side .....	40
Fig. 6: Division of the dimmable section into 200 brightness levels .....	41
Table 1: Pin assignment of the supply input socket .....	33
Table 2: Pin assignment of the "Trigger in" (input) .....	35
Table 3: Pin assignment of the "Sync out" (output) .....	35
Table 4: Pin assignment of the remote socket control voltage .....	35

# 1 FUNDAMENTALS

## 1.1 GENERAL INSTRUCTIONS


Familiarize yourself with the appliance and how to operate it by reading these operating instructions, before connecting the appliance with the mains and using it. The appliance must only be operated by duly instructed persons. Apart from that, the user must be familiar with the operation of the components connected with this appliance and how to connect these components. These operating instructions are part of the appliance as such and shall be retained over the appliance's entire service life. Should the user of the appliance change, these operating instructions must be passed on to the new user.

The appliance has been developed, manufactured and tested in accordance with EN standard 61010-1:2001 (IEC 1010-1) "General Safety Requirements for Electrical Equipment for Measurement, Control & Laboratory Use", with the relevant requirements of Annex 1 of Directive 2006/95/EC being observed.

The appliance satisfies all requirements set forth in Directive 2004/108/EC. The relevant standards observed in this appliance are listed in Section 3 "Technical data" of these operating instructions. The conformity with the standards referred to herein is also documented by the CE sign. These operating instructions use the following symbols to mark instructions and warnings:

**PLEASE NOTE**

This symbol points to an instruction which must be given special attention

**PLEASE NOTE**

De-energize the appliance, before carrying out the specified activity (unplug the appliance)

**PLEASE NOTE**

This symbol points to a danger that the appliance, or components connected with it, may be exposed to.

**PLEASE NOTE**

This symbol points to a danger that the user may be exposed to.

**PLEASE NOTE**

Danger caused by an electric shock.

## Fundamentals

General instructions || Instructions concerning the appliance's safety

LH-LED 100

### WARNING

Danger from high-energy light radiation. Do not look directly into the light beam. Wear goggles!

### WARNING

Danger from hot parts, wait until they have cooled down.

## 1.2 INSTRUCTIONS CONCERNING THE APPLIANCE'S SAFETY

The appliance must only be used for its intended purpose as described in these operating instructions. The manufacturer cannot be held liable, if the appliance is used for any other purposes not specified herein. In such case, all guaranty and warranty claims shall be forfeited. The manufacturer's liability as well as any guaranty and warranty claims shall also be forfeited, if damaged components, or components unsuitable for this appliance, are used, or if unauthorized service and repair works are carried out on the appliance.

If safety defects or malfunctions are identified, the appliance shall immediately be put out of operation and secured against any further use, with all poles being disconnected from the mains. The appliance shall be sent to the supplier, or any representative duly authorized by him, to have it repaired. The same shall apply to the power supply unit belonging to the appliance.

The appliance will generate high-energy light with invisible radiation. There is a danger of getting blinded and going blind. Never look straight into the appliance's light output. Always connect the appliance with the system to be lighted, before switching it on.

The power supply unit supplied with the appliance must only be plugged into a properly installed socket outlet with grounding contact. The power supply unit must only be operated in dry rooms and not in explosive ambiances.

Before connecting components with the appliance, disconnect all poles of the latter from the mains.

## Fundamentals

LH-LED 100

Instructions concerning the appliance's safety

The power supply unit supplied with the appliance is designed for a voltage range between 100 and 240 VAC and 50...60 Hz.

It will be automatically adjusted to the actual voltage within the above voltage and frequency range. When operating the appliance at a voltage outside this range, the appliance and the components connected with it may get destroyed.

## Fundamentals

Shipping, storage and unpacking instructions ||  
Warranty instructions

LH-LED 100

---

### 1.3 SHIPPING, STORAGE AND UNPACKING INSTRUCTIONS

The ambient conditions specified in these operating instructions for shipping and storing the appliance must be observed. The appliance shall be shipped in its original packaging only.

Should any damage be identified when unpacking the appliance, the supplier or the forwarder commissioned by him shall be immediately informed.

Check the completeness of all parts belonging to the light source's scope of delivery (see 2.3 "Scope of supply").

Let the appliance adjust itself to the ambient conditions for about 1 hour after the unpacking, before putting it into operation. This will prevent damage caused by condensate that may accumulate inside the appliance as a result of temperature differences.

### 1.4 WARRANTY INSTRUCTIONS

No other maintenance and repair works must be carried out on the appliance, except the ones listed in these operating instructions. Only the manufacturer of the appliance or duly authorized representatives are permitted to carry out repairs.

---

## 2 DESCRIPTION OF THE APPLIANCE

### 2.1 INTENDED USE

The lamp housing LH-LED 100 is a light source for being directly connected via manufacturer-specific adapters with microscopes from various manufacturers. Any other use is deemed to be not an intended purpose.

### 2.2 IDENTIFICATION OF THE APPLIANCE

The appliance's type plate with all details necessary to identify the appliance is affixed to its rear side.

### 2.3 SCOPE OF SUPPLY

The lamp housing's scope of supply includes the following components:

- lamp housing LH-LED 100
- power supply unit
- operating instructions
- collimator (optional, installed in the lamp housing)
- microscope adapter (optional)

---

### 2.4 SPECIAL FEATURES OF THE APPLIANCE

The LH-LED 100 can replace lamp housings that are operated with a 100-W halogen lamp. The LED's considerably higher service life as compared with that of a halogen lamp will provide the following benefits to the user:

- no change of lamps required
- no adjustment of lamp and reflector required
- higher availability of the entire system (microscope)
- no procurement of the replacement illuminant required

Besides, the lamp housing LH-LED 100 has only a maximum power input of 40 W as compared with a conventional 100-W halogen lamp housing.

The lamp housing's brightness can be adjusted at the appliance by pressing a button, on an operator console, via the external control voltage and via a serial interface (optionally via an adapter).

### 3 TECHNICAL DATA

Dimensions (width x height x depth) .....approx. 125 mm x 120 mm x 135 mm  
 Mass .....approx. 2 kg

**Ambient conditions**

Storage and shipment (packaged):

admissible ambient temperature.....-20 ... 85°C  
 admissible ambient temperature..... 10 to 90 %, non-condensing

Operation:

admissible ambient temperature..... 0 ... 40°C  
 admissible relative humidity ..... 10 to 90 %, non-condensing

**Operating data:**

operating ambience ..... dry and closed rooms  
 operating altitude ..... maximum 2 000 m a.s.l.  
 degree of protection ..... IP20 (EN 60529)

**EMC:**

.....DIN EN 61326-1 (2006-10)  
 This appliance satisfies the requirements of Directive 2004/108/EC.

**Electrical safety:**

Pursant to ..... EN 61010-1:2002

**Power supply:**

input voltage ..... 11,5 ... 12,5 V<sub>DC</sub>  
 power input ..... maximum 50 VA  
 supply input socket ..... 4-pole mini DIN Power Jack (such as Kycon KPJX-4)  
 Pin assignment ..... see Table 1

Connection	Function
1	+ U <sub>in</sub>
2	+ U <sub>in</sub>
3	GND
4	GND
shield	PE

Table 1 Pin assignment of the supply input socket

The lamp housing will be supplied with power from a separate power supply unit with the relevant isolation and approvals (included in the scope of supply). Using any other power supply units is not in conformity with the intended use.

Luminous flux @6.500 K maximum ..... 1.000 lm  
Compensation of the thermal drift of the luminous flux ..... < 5 %  
Color temperature ..... typically 6.500 K  
Spectral distribution ..... see diagram  
Average service life of the illuminant (luminous flux  $\leq 70\%$ ) ..... 50.000 h


Fig. 1 Spectral distribution of the luminous flux as compared with a halogen lamp

**Output (optical):**

Suitable for microscopes of the following manufacturers.. Zeiss, Leica, Nikon, Olympus  
via factory-specific adapters & collimators  
(optional accessories)

Adapters and collimators for other microscope manufacturers on request

**Interfaces:**

„Trigger in“ (input):

Function .....switching on and off the light source (LED)  
 Level ..... TTL (1 = off; 0 = on)  
 maximum trigger frequency ..... ≤ 25 kHz  
 Potential ..... potential-free (floating) related to the supply voltage  
 Socket.....BNC 50 Ohm  
 Pin assignment ..... see Table 2

Connection	Function
internal conductor	signal trigger
shield	signal trigger - GND

Table 2 Pin assignment of the “Trigger in” (input)

„Sync out“ (output):

Function .....signal synchronous PWM-Frequency  
 (H-level □ LED on)  
 Level ..... TTL  
 Potential ..... potential-free (floating) related to the supply voltage  
 Socket.....BNC 50 Ohm  
 Pin assignment ..... see Table 3

Connection	Function
internal conductor	signal out
shield	GND

Table 3 Pin assignment of the “Sync out” (output)

„Remote“ (input & output):

Socket..... 9pole D-Sub female  
 Function 1 ..... interface with the operator console  
 Function 2 ..... interface for the control voltage 0 ... 5V  
 Potential ..... related to the supply voltage GND  
 Pin assignment..... see Table 4

Connection	Function
pin 8	control voltage
pin 1	GND

Table 4 Pin assignment of the remote socket control voltage

Function 3 .....serial interface  
 (optional)

## 4 PUTTING THE APPLIANCE INTO OPERATION

### 4.1 POSITIONING THE APPLIANCE

Danger of blinding: The appliance will generate intensive light radiation in the visible spectral range. Always connect the lamp housing with the target system (such as a microscope) first.


Connect the lamp housing LH-LED 100 with the target system via a fitting adapter in such a way that the cooling fins are located in a vertical position and that the flow inside is not obstructed. Parts of the housing will warm up during operations and must not be covered.

Following that, the peripheral equipment and accessories (operator console, serial adapters...) can be connected. In a last step the lamp housing must be connected with the power supply unit and the latter with the mains. The cables connected with the lamp housing must not pose a danger.

### 4.2 FASTENING THE MICROSCOPE ADAPTER

The microscope adapter will be fastened by screwing it into the optical output of the lamp housing (Fig. 2). Before screwing in the adapter make sure that the collimator has been inserted into the lamp housing (a collimator is available separately). The lamp housing can now be attached to the microscope, but observe the instructions of the microscope manufacturer.

- 1 Optical output of the lamp housing
- 2 Microscope adapter (example)

Fig. 2: Fastening of the microscope adapter

### 4.3 CONNECTING ELECTRICAL EQUIPMENT

The lamp housing LH-LED 100 is equipped with several terminal connections, from which the lamp housing can be controlled (analog or digital control), or from which its status can be inquired. The lamp housing can thus be easily integrated into a system control.


Fig. 3: Electrical connections

- Remote socket (Fig. 3[1]): makes it possible to connect various accessories
  - Operator console (optional equipment): Allows the user to control the brightness of the lamp housing easily from his/ her workplace. In addition to that, the brightness currently set will be indicated on a display, while freely selectable brightness values can be chosen by pressing the relevant quick selection buttons.
  - Control voltage (analog): The brightness of the lamp housing can be adjusted via the control voltage that ranges between 0 and 5 V. When applying a voltage of 0 V, the light source of the lamp housing will be off, while it is shining with its full brightness at a voltage of 5 V. For the pin assignment for this function please refer to item "Interfaces-Remote-Function 2" in section 3 "Technical Data".
  - RS2RS 232 adapter/USB adapter (optional equipment): The brightness of the lamp housing can be adjusted via a control unit (such as a PC) with the help of these adapters.

## Commissioning

Connecting electrical equipment

LH-LED 100

- 
- “Sync out” (output) (Fig. 3[2]):  
The signal transmitted from this output will indicate the switching status of the light source in the lamp housing. This will make it possible, among other things, to synchronize the lamp housing with cameras, above all, when the lamp housing operates in the pulse width modulated status (approx.. 0 ... 30 % of brightness)  
For further details and the pin assignment for this function please refer to item “Interfaces - Sync out (output)” in section 3 “Technical Data”.
  - “Trigger in” (input) (Fig. 3[3]):  
The illuminant (i.e. the LED) of the lamp housing can be switched on and off via this input.  
For further details and the pin assignment for this function please refer to item “Interfaces – Trigger in (input)” in section 3 “Technical Data”.
  - Socket for the supply voltage (input) (Fig. 3[4]):  
Connection for the power supply unit that is included in the scope of supply. Using any other power supply unit is not in conformity with the intended use.

## 5 OPERATION

### 5.1 OPERATING AND DISPLAY ELEMENTS


Fig. 4: View at the lamp housing's operating & display elements

As soon as the lamp housing is supplied with power, the indicator in the ON/OFF button (Fig. 4[1]) will light up in red. On pressing this button, the lamp housing will be switched on and the indicator in the button will light in green.

After the lamp housing has been switched on, the LED (i.e. the light source) will shine with minimum brightness (glaring protection). The brightness of the LED can now be adjusted by pressing the "-" or "+" button (Fig. 4[2 & 3]).


Alternatively, the brightness of the LED can also be adjusted via the control voltage which is applied through the "Remote" socket (Fig. 3[1]). If the "-" or "+" button is now pressed, the brightness level set during such action will then be accepted as the default value and changes to the lamp housing can only be made by pressing the "-" or "+" buttons. The control voltage can only be used as a means of controlling the brightness again, after the ON/OFF button has been briefly pressed. The lamp housing will be automatically switched on, as soon as the control voltage is applied.

An operator console, which is optionally available, will make it possible to adjust the LED's brightness easily from the user's workplace. Besides, the operator console also allows the storage of freely selectable brightness values and their choice by pressing the relevant key just once.

As a further option, RS232 or USB adapter cables are also available that will then make it possible to adjust the LED's brightness from a PC.

### 5.2 COLLIMATOR ADJUSTMENT

The adjusting device for the collimator has been installed at the lamp housing's bottom side. It can be used to adapt the divergence of the light bundle to the target system (microscope). After loosening the fastening bolt (Fig. 5[1]) with a hexagon socket wrench SW 1.5 (Fig. 5[2]) by about half a turn, the collimator can be (re-)set by shifting the slide along the lamp housing's optical axis and re-arresting it.

*Fig. 5: Collimator adjustment on the bottom side*

### 5.3 SETTING THE BRIGHTNESS LEVEL

The dimming is achieved at a constant LED light flux up to about 30 % of the maximum brightness by pulse width modulation (PWM). The modulation curve for the synchronization of external (peripheral) equipment (such as cameras) will be provided as TTL signal from the “Sync out” output. Above 30 % of the brightness level, the light flux will be adjusted and set by the LED in 100 steps and no longer be modulated.


Fig. 6: Division of the dimmable section into 200 brightness levels

*Steps 0 ... 100 (corresponds to the control voltage of 0 ... 2.5 V):*

*Dimming through PWM at 23.9 kHz, constant LED light flux*

*Steps 100 ... 200 (corresponds to the control voltage of 2.5 ... 5 V):*


*Dimming through changing the LED light flux without PWM*

## 6 SERVICE

### 6.1 MAINTENANCE AND CLEANING

The lamp housing LH-LED 100 does not require regular maintenance. This does not release the user from his duty to have it inspected regularly in accordance with the statutory regulations for the purpose of accident protection and to ensure the appliance's electrical safety.

Before attempting any cleaning work, switch off the lamp housing and let it cool down. Use a dry, lint-free cloth only. When cleaning the appliance make sure that no dust or other foreign particles can get inside through the ventilation openings. Should the lamp housing be damaged in any way, also externally, take it out of operation immediately and return it for servicing.

Should the lamp housing LH-LED 100 not be used for a longer period, protect it with a suitable cover against dust (e.g. the plastic bag of the original packaging).

### 6.2 SPARE PARTS & ACCESSORIES

The lamp housing LH-LED 100 has not any wear and spare parts that the user might be required to exchange. Only the optional equipment, such as the microscope adapter and the collector, can be additionally fitted and replaced, if the need arises.

### 6.3 REPAIRS

Repairs must only be carried by the manufacturer of the appliance or a duly authorized representative

Should the lamp housing LH-LED 100 be defective, so that it needs to be returned to the supplier or his duly authorized representative, use the original packaging for the shipment.


Adapter Nikon Bayonet  
57 05 01-001-27

Adapter Nikon  
55 01 02-804-10



Adapter Olympus  
55 01 02-803-10


Dovetail LED-Z  
55 01 02-802-10

Dovetail LED-L  
55 01 04-801-10


Operator Console LH-LED  
57 04 01-001-26


USB-Adapter  
57 04 11-001-27


Lamp Housing LH-LED 100  
55 01 02-001-26

Collimator  
55 01 02-200-14

Mains Connecting Cable EU  
(Plug-Type E+F)

Mains Connecting Cable  
other Regions


Table Power Unifitil  
(Included in scope of delivery)